


# TQS Electronic Tourniquet

## Specification:

- Quick, Quiet and easy to use
- Single Clear LCD, Dual channel Pre-set and cuff pressure Display
- Designated support LCD with controls
- SURTRAK fully, easy programable audio and visual surgical time tracking
- Programable Quick Pressure pre-set function Button
- On Board Memory, surgical procedure recording and PC software
- Utility cart or drip pole mountable
- Audio and visual alarms



**dispozee**

## Disposable Tourniquet Cuffs

Tailor your sample cuffs with Oak Medical Services Dispozee Cuff Range:



**Fastener?**

- Velcro
- Adhesive Tape

**Ties?**

- Velcro
- Ties

**Fit?**

- Straight
- Conical

**Style?**

- Disposable
- 10 Use
- Disposable Sleeve

**dispozee**  
**Prevee-prep**  
for use with disposable & re-useable tourniquet cuffs



**OAK MEDICAL**  
SERVICES LTD

Oak Medical Services Ltd  
Unit 5A Albert Street  
Brigg, North Lincolnshire DN20 8HQ

T: 01652 657200 F: 01652 657009

W: [www.oakmedicalservices.co.uk](http://www.oakmedicalservices.co.uk)

E: [info@oakmedicalservices.co.uk](mailto:info@oakmedicalservices.co.uk)

## Power-Pod

Try Oak Medical Services range of "Power-Pod" units!

- 6-10 gang sockets version available
- Lockable castors
- Service and testing

